

BIO PLAS, Inc.

Pestles

BIO PLAS disposable polypropylene pestles may be used to homogenize cells or re-suspend precipitates such as protein or DNA pellets in buffer solutions. All the packaging options contain products free of detectable RNase, DNase, DNA and Pyrogens. These pestles and tubes are precision manufactured and quality controlled to firmly match together for maximum performance efficiency. Pestles are conveniently packaged into the appropriate size tubes, ready for use. All pestles are RNase/DNase free and sterilized.



*BIO PLAS Pestle
with 1.5mL "G"-Tube™*



*BIO PLAS Pestle
with 0.6mL "G"-Tube™*

Ordering Information

BPI Cat #	Tube Size	Pestle Color	Qty/Pack	Packs/Case
4030-PB	1.5mL	Blue	10 per Bag	25 Bags
4040-PB	0.6mL	Blue	10 per Bag	25 Bags

NOTE: *Pestles and G-Tubes packaged with pestles are a non Siliconized product*

Minimum order one case.