

BIO PLAS Inc.

Histology Specimen Protection

Bio Plas, Inc. is pleased to offer you both Foam Biopsy Pads and our new Nylon Biopsy Bags for all your Histology Specimen Protection!

Foam Biopsy Pads

The FOAM BIOPSY PADS are fabricated from cellular polyester plastic foam which is manufactured uniformly with air pockets throughout to allow for adequate fluid movement around the tissue while it is securely sandwiched between two pads. This assures the orientation of small biopsy tissue. Available in sizes to fit TISSUE PROCESSING CAPSULES and CAPSETTES™



Nylon Biopsy Bags

Our Nylon Biopsy Bags protect small specimens during processing and allow full penetration of reagents around tissue, and reduce the reagent carry over. Once processing is complete, bags provide quick access by simply peeling open for retrieval.

ORDERING INFORMATION

Catalog #	Color	Description	Qty/Pack	Packs/Case
6100	Blue	Foam Biopsy Pads (Use with 38mm capsules)	1,000	10
6110	Blue	Foam Biopsy Pads (Use with 29mm capsules)	1,000	10
6120	Blue	Foam Biopsy Pads (Use with 1 x 1-¼ cassettes)	1,000	10
6130	White	Foam Biopsy Pads (Use with 1 x 1-¼ cassettes)	1,000	10
6200	White	Nylon Biopsy Bags - 1" x 2"	100	10
6222	White	Nylon Biopsy Bags - 1.22" x 1.65"	100	10