


BIO PLAS

Bacti Cell Spreaders 30mm & 60mm Wide

- Autoclavable, non-breakable (eliminates glass hazard)
- 30mm head for prepared bi-plates
- 60mm head for 100-150mm prepared plates
- Sturdy plastic provides uniform design
- Easy to manipulate without slippage
- Easy to manipulate around curves and into corners of prepared plates
- Avoids waste of “disposables” thus becoming an environmental benefit
- Always prevents agar surface tears
- Allows for even distribution of specimen
- Produces an even smear for any permanent slide staining, i.e., PVA preserved fecal specimen for protozoan parasites
- Microbiological research of all types, especially for distribution of dilutions, i.e., colony counts.



The BIO PLAS BACTI CELL SPREADERS are unique products designed to greatly simplify and speed up bacterial culture work.

There are two sizes: 30mm and 60mm wide, which allow these polypropylene BACTI CELL SPREADERS to make it very easy to evenly distribute bacterial cells rapidly over an entire nutrient surface.

The BIO PLAS BACTI CELL SPREADERS are individually wrapped (or available in bulk) and may be cleaned and autoclaved for multiple uses.

Ordering Information

Catalog Number	Description	Pack Qty.	Packs/Case
6400	30mm BACTI CELL SPREADER, Blue	25	4
6401S	30mm BACTI CELL SPREADER, Blue, Sterilized	25 Individually Wrapped	4
6405	60mm BACTI CELL SPREADER, Yellow	25	4
6406S	60mm BACTI CELL SPREADER, Yellow, Sterilized	25 Individually Wrapped	4

BIO PLAS, inc., 4340 Redwood Highway, Suite A1, San Rafael, CA 94903 USA

(415) 472-3777 – Tel (415) 472-3758 – Fax www.bioplas.com – Web Site info@bioplas.com - Email