


BIO PLAS, Inc.

Polyester Foam Test Tube & Microcentrifuge Tube Racks

Features:

- Flame Retardant
- Color – Charcoal
- Long wearing
- Chemical and solvent resistant
- Float on H₂O bath with full tube and sample
- Ideal for refrigerated or frozen storage samples
- All tubes can be decanted simultaneously
- Low profile- will not tip over easily
- Lightweight- weighs less than two ounces
- Holds and locks 50 culture tubes securely
- Holds and locks 30 microcentrifuge and micro blood collecting tubes securely in racks
- More uniform incubation- tubes may be pushed down into H₂O bath while still being held securely in rack
- May be used in blood drawing trays-holds blood drawing tubes securely



Ordering Information

BPI Cat. #	Tube Size	Rack Capacity	Pkg. Qty.	Pkgs/Case
0010	10 x 75mm	50 Tubes	6 racks	5
	12 x 75mm	50 Tubes	6 racks	5
	13 x 100mm	50 Tubes	6 racks	5
0011	15 x 85mm	50 Tubes	6 racks	3
	16 x 100mm	50 Tubes	6 racks	3
	16 x 125mm	50 Tubes	6 racks	3
	16 x 150mm	50 Tubes	6 racks	3
0020	250µl	30 Tubes	6 racks	6
	400µl	30 Tubes	6 racks	6
0021	500µl	30 Tubes	6 racks	6
	1.5ml	30 Tubes	6 racks	6
	1.8ml	30 Tubes	6 racks	6