

Bio Plas, Inc.

Micro Tube Tool

This handy polypropylene tool opens Snap-Cap and Screw-Cap Microcentrifuge Tubes easily, prevents fatigue and saves wear on your fingers and thumbs!



The Micro tube tool's larger end is a socket wrench for Screw-Cap Micro Tubes.



Shown: The Micro Tube tool's smaller end is used to open a Snap Cap Micro Tube.

The innovative Micro Tube Tool makes it a breeze to open tight-fitting Snap Caps -and it lets you remove Screw Caps while tubes are still in the rack. Convenient and inexpensive, the Micro Tube Tool puts an end to difficult maneuvering and straining to open Screw-Cap Tubes and any size Snap Cap Tubes.

ORDERING INFORMATION

Cat #	Description	Pkg. Qty.	Pkgs./ Case
4000	Micro Tube Tool	4	25

BIO PLAS, INC

4340 Redwood Highway, Suite A1, San Rafael, CA 94903, (415) 472-3777 Tel / (415) 472-3758 Fax

www.bioplas.com – info@bioplas.com