


Snap-Cap Microcentrifuge Tubes



Solid sealing in a super-tough tube

Five sizes: 250 μ l, 400 μ l, 500 μ l, 1.5ml and 1.8ml. All sizes have attached, tight-fitting Snap-Caps for optimal sealing during sample processing and centrifugation. Colors: Natural (white) for 250 μ l, 400 μ l, and 500 μ l. The 1.5 ml and 1.8ml sizes are available in natural and four additional colors for color-coding.

Our Snap-Cap tube gives a firm, dependable seal, and its tough virgin-polypropylene construction withstands the rugged G-forces of high-speed centrifugation. Use it with confidence for specimen collection, mixing, electrophoresis sampling/handling, enzyme assays, radioassays, immunoassays and many other applications.

The Snap-Cap tube's smooth, conical interior lets you remove your entire sample with ease, and without material hangup. You can even cut it to remove a spun pellet. The Snap-Cap is attached to each tube, and the entire unit is autoclavable at 250°F at 15 psi for 15 minutes. Fits Beckman, Eppendorf and all other microcentrifuges.

Ordering Information

Cat #	Description	Color	Pkg. Qty.	Pkgs./Case
4003	Micro-Centrifuge Tube 400 μ l, Polypropylene	White	1000	10
4004	Micro-Centrifuge Tube 400 μ l, Polyethylene	White	1000	10
4005	Micro-Centrifuge Tube 500 μ l, Polypropylene	White	1000	10
4025	Micro-Centrifuge Tube 250 μ l, Polyethylene	White	1000	10
4150	Micro-Centrifuge Tube 1.5ml, Polypropylene	White	1000	10
4151	Micro-Centrifuge Tube 1.5ml, Polypropylene	Red	1000	10
4152	Micro-Centrifuge Tube 1.5ml, Polypropylene	Yellow	1000	10
4153	Micro-Centrifuge Tube 1.5ml, Polypropylene	Blue	1000	10
4154	Micro-Centrifuge Tube 1.5ml, Polypropylene	Green	1000	10
4155	Micro-Centrifuge Tube 1.8ml, Polypropylene	White	1000	10
4156	Micro-Centrifuge Tube 1.8ml, Polypropylene	Red	1000	10
4157	Micro-Centrifuge Tube 1.8ml, Polypropylene	Yellow	1000	10
4158	Micro-Centrifuge Tube 1.8ml, Polypropylene	Blue	1000	10
4159	Micro-Centrifuge Tube 1.8ml, Polypropylene	Green	1000	10

Micro Tube Tool

This handy polypropylene tool opens Snap-Cap and Screw-Cap Microcentrifuge Tubes easily, prevents fatigue and saves wear on your fingers and thumbs!

The innovative Micro Tube Tool makes it a breeze to Open tight-fitting Snap Caps- and it lets you remove Screw Caps while tubes are still in the rack. Convenient and Inexpensive, the Micro Tube Tool puts an end to difficult Maneuvering and straining to open Screw-Cap Tubes and any size Snap-Cap Tube.



Ordering Information

Cat #	Description	Pkg. Qty.	Pkgs./ Case
4000	Micro Tube Tool	4	25



BIO PLAS, INC

4340 Redwood Highway, Suite A1
San Rafael, CA 94903
(415) 472-3777 Tel / (415) 472-3758 Fax
www.bioplas.com – info@bioplas.com