

Uni-Flex™ Safety Caps

Now there's a safe, economical way to cap and recap blood collection tubes and disposable glass culture tubes, regardless of manufacturer. Uni-Flex™ Safety Caps are available for a wide range of biomedical and industrial applications.

The Uni-Flex™ Safety Cap liquid tight seal ensures sample integrity and prevents spillage of serum samples. Uni-Flex™ Safety Caps are available in several colors which helps guard against sample cross-contamination and is convenient for color-coding to particular R&D projects. Uni-Flex™ Safety Caps protect your laboratory staff from the aerosols of highly infectious organisms such as TB and the HTLV-III virus.

- Will not aerosol when opened!¹
- 10, 12, 13 and 16mm sizes fit all varieties of blood collecting and disposable culture tubes.
- Liquid tight seal prevents sample evaporation. The Uni-Flex™ Safety Cap is designed to hold under centrifugation, agitation, refrigeration and freezing.
- Easy identification of samples. Uni-Flex™ Safety Caps are available in seven colors that conveniently correspond to blood collection cap colors.
- Stretchable, unique over-locking design is far superior to Parafilm™. Locks securely onto the outside of the tube lip to prevent evaporation.
- Easy, one hand removal and application. To close, press down firmly with thumb to lock in place. To open, simply pry cap up with thumb.
- Easy to label. Write directly on the cap with indelible ink.
- Clear, natural colored caps are trace-metal free.
- OEM product packaging.



Using a Uni-Flex™ Safety Cap is an easy, one hand operation. To close, press down firmly with thumb to lock in place.

To open, simply pry cap up with thumb.



References

¹ "HIV Infection", MMWR, Vol. 36/No. 19, May 22, 1987

"What every laboratorian should know about AIDS", Medical Laboratory Observer, Nov., 1985.

"The AIDS Virus", New York Times, May 2, 1987

Uni-Flex™ Safety Cap is a trademark of BIO PLAS, Inc., San Rafael, CA

Parafilm is a trademark of American Can Co./Dixie Marathon.



BIO PLAS, Inc.

4340 Redwood Highway, Suite A1
San Rafael, CA 94903 USA

(415) 472-3777 – Tel / (415) 472-3758 – Fax

www.bioplas.com

info@bioplas.com