


## Polyester Foam Test Tube & Microcentrifuge Tube Racks

### Features:

- Flame Retardant
- Color – Charcoal
- Long wearing
- Chemical and solvent resistant
- Float on H<sub>2</sub>O bath with full tube and sample
- Ideal for refrigerated or frozen storage samples
- All tubes can be decanted simultaneously
- Low profile- will not tip over easily
- Lightweight- weighs less than two ounces
- Holds and locks 50 culture tubes securely
- Holds and locks 30 microcentrifuge and micro blood collecting tubes securely in racks
- More uniform incubation- tubes may be pushed down into H<sub>2</sub>O bath while still being held securely in rack
- May be used in blood drawing trays-holds blood drawing tubes securely



### Ordering Information

| BPI Cat. # | Tube Size  | Rack Capacity | Pkg. Qty. | Pkgs/Case |
|------------|------------|---------------|-----------|-----------|
| 0010       | 10 x 75mm  | 50 Tubes      | 6 racks   | 5         |
|            | 12 x 75mm  | 50 Tubes      | 6 racks   | 5         |
|            | 13 x 100mm | 50 Tubes      | 6 racks   | 5         |
| 0011       | 15 x 85mm  | 50 Tubes      | 6 racks   | 3         |
|            | 16 x 100mm | 50 Tubes      | 6 racks   | 3         |
|            | 16 x 125mm | 50 Tubes      | 6 racks   | 3         |
|            | 16 x 150mm | 50 Tubes      | 6 racks   | 3         |
| 0020       | 250µl      | 30 Tubes      | 6 racks   | 6         |
|            | 400µl      | 30 Tubes      | 6 racks   | 6         |
| 0021       | 500µl      | 30 Tubes      | 6 racks   | 6         |
|            | 1.5ml      | 30 Tubes      | 6 racks   | 6         |
|            | 1.8ml      | 30 Tubes      | 6 racks   | 6         |

# BIO PLAS, INC.


## 96 and 64 Well "Reversible" Micro Tube Racks and Snap-On Rack Covers



Both the 96 and 64 racks are designed so that one side of the rack accommodates 0.5ml microcentrifuge tubes. The reverse side holds 1.5ml to 2.0ml microcentrifuge tubes and screw cap microcentrifuge tubes.

The clear polypropylene cover on our product easily snaps on to keep it firmly in place during transportation and storage. A recessed area in the top of the covers allows for secure stacking of racks one on top of each other. The snap on cover is designed to accommodate the extra height of our screw cap microcentrifuge tubes.

- Available in a wide variety of colors
- Two frosted writing areas for sample identification
- Autoclavable: 15 minutes, 15lbs. @ 120° C.
- Will withstand freezer storage to -90°C
- Alpha-numeric markings, "A" thru "H" vertically, 1 thru 8 horizontally
- Replaces the 5-1/4" x 5-1/4" freezer storage box



To the Left: BIO PLAS INC. 64 Well "Reversible" Micro Tube Racks, Above: BIO PLAS INC. 96 Well "Reversible" Micro Tube Racks.

### Ordering Information

| 64 Well<br>Cat # | 96 Well<br>Cat # | Description     | Color                 | Pkg.<br>Qty. | Pkgs.<br>Case |
|------------------|------------------|-----------------|-----------------------|--------------|---------------|
| 0080             | 0090             | Micro Tube Rack | Natural               | 5            | 4             |
| 0081             | 0091             | Micro Tube Rack | Assorted*             | 5            | 4             |
| 0082             | 0092             | Micro Tube Rack | Blue                  | 5            | 4             |
| 0083             | 0093             | Micro Tube Rack | Green                 | 5            | 4             |
| 0084             | 0094             | Micro Tube Rack | Orange                | 5            | 4             |
| 0085             | 0095             | Micro Tube Rack | Pink                  | 5            | 4             |
| 0086             |                  | Micro Tube Rack | Yellow                | 5            | 4             |
| 0087             |                  | Micro Tube Rack | Lavender              | 5            | 4             |
| 0089             | 0099             | Micro Tube Rack | Black                 | 5            | 4             |
| 0081F            | 0091F            | Micro Tube Rack | Fluorescent Assorted* | 5            | 4             |
| 0082F            | 0092F            | Micro Tube Rack | Fluorescent Blue      | 5            | 4             |
| 0083F            | 0093F            | Micro Tube Rack | Fluorescent Green     | 5            | 4             |
| 0084F            | 0094F            | Micro Tube Rack | Fluorescent Orange    | 5            | 4             |
| 0085F            | 0095F            | Micro Tube Rack | Fluorescent Pink      | 5            | 4             |
| 0086F            | 0096F            | Micro Tube Rack | Fluorescent Yellow    | 5            | 4             |
| 0050             | 0052             | Tube Rack Cover | Natural               | 5            | 4             |


\* Assorted packages consist of one each: Blue, Green, Orange, Pink and Yellow

Min. Order One Package

# BIO PLAS, INC.

## 80 Well Micro Tube Racks and Rack Covers

The 80 Well Micro Tube Racks are designed with 80 wells in a 16 x 5 configuration with each well numbered and two textured writing surfaces allowing convenient sample identification. Rack dimensions are 9"L x 2 1/16"D with a well diameter of 10.5mm.



Made of polypropylene, the 80 Well Micro Tube Racks are autoclavable, ideal for freezing, benchtop use, storage and transporting samples.

Available in 9 original colors plus 5 new fluorescent colors, this assortment of colors allows lab technologists to color code specific projects and allows for identification of hazardous experiments.

**May be used conveniently with the following size microcentrifuge tubes:**

| <u>Snap Cap MCT's</u> | <u>Screw Cap MCT's</u> |
|-----------------------|------------------------|
| ◦ 500ul               | ◦ 0.5ml                |
| ◦ 1.5ml               | ◦ 1.5ml                |
| ◦ 1.8ml               | ◦ 2.0ml                |

The rack covers are manufactured with clear polypropylene and easily snap onto the 80 Well Micro Tube Racks. Cover tops have an indented area that allows for secure stacking. Both racks and covers are autoclavable at 252° (121°C) for 15 minutes at 15 pounds.

## Ordering Information

| BPI Cat # | Description             | Color                 | Pkg. Qty. | Pkgs/Case |
|-----------|-------------------------|-----------------------|-----------|-----------|
| 0060      | 80 Well Micro Tube Rack | Natural               | 5         | 4         |
| 0061      | 80 Well Micro Tube Rack | Assorted*             | 5         | 4         |
| 0062      | 80 Well Micro Tube Rack | Blue                  | 5         | 4         |
| 0063      | 80 Well Micro Tube Rack | Green                 | 5         | 4         |
| 0064      | 80 Well Micro Tube Rack | Orange                | 5         | 4         |
| 0065      | 80 Well Micro Tube Rack | Pink                  | 5         | 4         |
| 0066      | 80 Well Micro Tube Rack | Yellow                | 5         | 4         |
| 0067      | 80 Well Micro Tube Rack | Lavender              | 5         | 4         |
| 0068      | 80 Well Micro Tube Rack | Red                   | 5         | 4         |
| 0069      | 80 Well Micro Tube Rack | Black                 | 5         | 4         |
| 0070      | 80 Well Micro Tube Rack | Fluorescent Blue      | 5         | 4         |
| 0071      | 80 Well Micro Tube Rack | Fluorescent Green     | 5         | 4         |
| 0072      | 80 Well Micro Tube Rack | Fluorescent Orange    | 5         | 4         |
| 0073      | 80 Well Micro Tube Rack | Fluorescent Pink      | 5         | 4         |
| 0074      | 80 Well Micro Tube Rack | Fluorescent Yellow    | 5         | 4         |
| 0075      | 80 Well Micro Tube Rack | Fluorescent Assorted* | 5         | 4         |
| 0051      | Tube Rack Cover         | Natural               | 5         | 4         |

\* Assorted packages consist of one each: Blue, Green, Orange, Pink and Yellow.

Minimum Order One Package



### BIO PLAS, INC

4340 Redwood Highway  
Suite A15

San Rafael, CA 94903

(415) 472-3777 Tel / (415) 472-3758 Fax

[www.bioplas.com](http://www.bioplas.com) – [info@bioplas.com](mailto:info@bioplas.com)