



Screw Cap Microcentrifuge Tubes Material Safety Data Sheet

Product Name: Screw Cap Microcentrifuge Tubes

Bio Plas Catalog Numbers : 4200 through 4224 & 4200SB through 4204ASB

Company Name: Bio Plas, Inc.

Address: 4380 Redwood Highway, Suite C-16, San Rafael, CA 94903

Tel: (415) 472-3777

Recommended Use: Sample Collection And Microcentrifugation

Product Composition: Polypropylene. Free of RNase, DNase, DNA, Pyrogens and PCR inhibitors.

Hazards Identification: Not classified as dangerous for use/supply

Hazardous Ingredients: None

Fire Hazard: Product burns at 390 Degrees Celsius. The fumes are CO and CO2

Fire Extinguishing Media: Water, Foam, Carbon Dioxide, Dry Chemical

Toxicology Information: Material is non-toxic

Storage Requirements:

- This product may react with strong oxidising agents & should not be stored near such materials
- Storage temperature should be ambient. (preferably below 50 degrees C)
- Open flames prohibited.
- Protect from heat and direct sunlight.

Note: Bio Plas Screw Cap Tubes are only guaranteed to work with Bio Plas Screw Caps.