



Bio Plas Capsettes & Uni-Capsettes Material Safety Data Sheet

Product Name: Capsettes And Uni-Capsettes

Bio Plas Catalog Numbers: 6030 through 6045 & 6050 through 6065

Company Name: Bio Plas, Inc.

Address: 4380 Redwood Highway, Suite C-16, San Rafael, CA 94903

Tel: (415) 472-3777

Recommended Use: Processing, embedding, sectioning, and storage of specimen samples.

Product Composition: Polypropylene. Free of RNase, DNase, DNA, Pyrogens and PCR inhibitors.

Hazards Identification: Not classified as dangerous for use/supply

Hazardous Ingredients: None

Fire Hazard: Product burns at 390 Degrees Celsius. The fumes are CO and CO2

Fire Extinguishing Media: Water, Foam, Carbon Dioxide, Dry Chemical

Toxicology Information: Material is non-toxic

Storage Requirements:

- This product may react with strong oxidizing agents & should not be stored near such materials
- Storage temperature should be ambient. (preferably below 50 degrees C)
- Open flames prohibited.