



Bio Plas, Inc.

New Improvement

Cat# 3250 - 5 ML MACRO TIPS

Product Change Notice!

SUMMARY OF CHANGE


- We have incorporated a new mold for these tips to improve our efficiencies and increase our production capacity. There is a slight change in the number of ribs on the top exterior crown of the tip which has no effect on its fit, form or function.

WHY ARE WE MAKING THIS CHANGE?

- The new mold incorporates a hot-runner system that facilitates faster cycle times and reduces plastic waste by 25%, compared to the cold-runner system.

HOW WILL THIS CHANGE AFFECT MY CUSTOMERS?

- Your customers may notice the change in the appearance of the 3250 tips. Again, only the ribbing is changed; the tips are unchanged in fit, form and function.



Dimensions are in millimeters.

Cat#	Change Type	Description	Packaging
3250	Cosmetic / Ribbing	5 mL Macro Tip, Eppendorf-Fit, Bulk	250 tips/unit, 10 units/case

If you have questions on this notification, please do not hesitate to contact us at info@bioplas.com.

Regards,

Bill Slattery
Bio Plas, Inc.