

## **Individual PCR Tubes**

- Individual PCR tubes feature evaporation-resistant, easy-to-open-and-close lids that prevent sample loss.
- Ultra-thin, uniform wells ensure optimal heat transfer and reaction efficiency.
- Conveniently packaged in reclosable bags to minimize risk of contamination.
- Available in ultra-clear and assorted.
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.

## PCR Strip Tubes And Strip Caps

- Robust linkages make the strips rigid and easy to handle, preventing spills.
- Angularly attached caps prevent lids and hinges from interfering with each other, allowing for easier operation in an SBS format rack.
- Unique two-material white tubes with optically-clear flat caps create the perfect environment for consistent, reliable qPCR performance. White wells offer enhanced fluorescence reflection and optically-clear caps ensure that the fluorescence signal is uninhibited.
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.





Angularly attached caps on BioFlex 8-Strip and Snappable tubes prevent cap lids and hinges from interfering with each other, allowing for easier operation in an SBS format rack.



Bio Plas, Inc. A Woman Owned Small Business www.bioplas.com info@bioplas.com Tel: (415) 472-3777 Fax: (415) 472-3758





## Flat Cap Hinged PCR Tubes & Caps

- Strips of hinged attached caps prevent accidents common with separate caps.
- Optically-clear, flat caps for qPCR.
- Uniform thin walls.
- For use in pharmaceutical industry, food industry, and molecular biology.
- Manufactured from high-quality, virgin polypropylene to ensure purity and provide excellent optical clarity.
- Certified free of detectable RNase, DNase, DNA, PCR inhibitors, and tested pyrogen-free.



Cat# 5060-10



Cat# 5060-9

Ordering Information				
Catalog #	Description	Color	Pkg. Qty	Pkgs./ Case
5060-1	BioFlex Individual PCR Tube with Frosted Flat Cap, 0.2ml	Clear	1000	10
5060-2	BioFlex Individual PCR Tube with Frosted Flat Cap, 0.2ml	Asst	1000	10
5060-3	BioFlex 8-Strip Tubes with Individually-Attached, Angular, Flat Caps, 0.2ml	Clear	120	10
5060-4	BioFlex 8-Strip Tubes with Individually-Attached, Angular, Flat Caps, 0.2ml	Asst	120	10
5060-5	BioFlex 8-Strip Low-Profile* Tubes with Individually-Attached, Angular, Flat Caps, 0.1ml	Clear	120	10
5060-6	BioFlex 8-Strip Low-Profile* Tubes with Individually-Attached Flat, Angular, Caps, 0.1ml	White	120	10
5060-7	BioFlex 8-Strip PCR Tubes & 8-Strip Flat Caps, 0.2ml	Clear	125	10
5060-8	BioFlex 8-Strip Low-Profile PCR Tubes & 8-Strip Flat Caps, 0.1ml	Clear	125	10
5060-9	BioFlex Snappable 8-Strip Tubes with Individually-Attached Flat Caps, 0.2ml	Clear	120	10
5060-10	BioFlex 8-Strip Tubes with Hinged Attached Flat Caps, 0.2ml	Clear	125	10

\* Fast<sup>®</sup> Cycler and Lightcycler<sup>®</sup> compatible.

Bio Plas, Inc. A Woman Owned Small Business www.bioplas.com info@bioplas.com Tel: (415) 472-3777 Fax: (415) 472-3758