



BIO PLAS INC.

Premium Laboratory Consumables
Made In The USA



Bio Plas Micro Tube Tool

This handy polypropylene **Micro Tube Tool** opens Snap-Cap and Screw Cap Microcentrifuge Tubes easily, prevents fatigue and saves wear on your fingers and thumbs!

The innovative **Micro Tube Tool** makes it a breeze to open tight-fitting Snap Caps- and it lets you remove Screw Caps while tubes are still in the rack. Convenient and inexpensive, the **Micro Tube Tool** puts an end to difficult maneuvering and straining to open Screw-Cap Tubes and any size Snap Cap Tubes.



Ordering Information			
Catalog #	Description	Pkg. Qty	Pkgs./Case
4000	Micro Tube Tool	4	25

Micro Tube Tool, Patent # 5,253,551



The **Micro Tube Tool's** larger end is a socket wrench for Screw-Cap Micro Tubes.

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA

Tel: (415) 472-3777

Fax: (415) 472-3758

info@bioplas.com

www.bioplas.com