



# BIO PLAS INC.

Premium Laboratory Consumables  
Made In The USA



## Histo Plas: Capsettes™ & Uni Capsettes™

Histo Plas™ Capsettes™ and Uni-Capsettes™ provide a modular, organized means for safe and consistent processing, embedding, sectioning and storage of specimen samples. Available in 16 colors which provide a convenient means of separating surgical and autopsy cases, institutions, pathologist, etc. Both the Histo Plas Capsettes and Uni-Capsettes are accepted by all of the current automatic cassette labelers.

Cost effective, yet quality engineered to remain sealed during processing. Strategically designed slots allow for efficient solvent flow and drainage, Histo Plas Capsettes and Uni-Capsettes will not float during processing. Our rigid design eliminates undersized cassettes during sectioning. Proprietary acetyl copolymer plastic is utilized to resist heat, solvents, decalcifying solutions, and problems associated with microwave processing.

Histo Plas Capsettes and Uni-Capsettes will fit standard stainless steel, and disposable plastic base molds. Uni-Capsettes are sold in packages of 500 disposable plastic tops and (500 specimen) bottoms. Technologists no longer need to struggle with detachment of the lid during embedding.

**Available in 16 Colors!**



Capsettes™, Blue, 6034



Uni-Capsettes™, Aqua, 6059



Capsettes™, Pink, 6032

**BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS**  
4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA  
Tel: (415) 472-3777 Fax: (415) 472-3758  
info@bioplas.com www.bioplas.com

Histo Plas Capsettes™ and Uni-Capsettes™ are trademarks of Bio Plas, Inc., San Rafael, CA



**Histo Plas™ Capsettes™ and Uni-Capsettes™ Benefits:**

- Heavy wall Construction will not float
- Will not open during processing
- Economical, 16 Colors
- Engineered to eliminate “chatter” during sectioning



Uni-Capsettes™, Burnt Orange 6056

Ordering Information				
Uni-Capsettes™ Catalog #	Capsettes™ Catalog #	Description	Pkg. Qty	Pkgs./ Case
6050	6030	Uni-Capsettes™/Capsettes™, White	500	10
6051	6031	Uni-Capsettes™/Capsettes™, Yellow	500	10
6052	6032	Uni-Capsettes™/Capsettes™, Pink	500	10
6053	6033	Uni-Capsettes™/Capsettes™, Green	500	10
6054	6034	Uni-Capsettes™/Capsettes™, Blue	500	10
6055	6035	Uni-Capsettes™/Capsettes™, Tan	500	10
6056	6036	Uni-Capsettes™/Capsettes™, Burnt Orange	500	10
6057	6037	Uni-Capsettes™/Capsettes™, Lilac	500	10
6058	6038	Uni-Capsettes™/Capsettes™, Gold	500	10
6059	6039	Uni-Capsettes™/Capsettes™, Aqua	500	10
6060	6040	Uni-Capsettes™/Capsettes™, Red	500	10
6061	6041	Uni-Capsettes™/Capsettes™, Chartreuse	500	10
6062	6042	Uni-Capsettes™/Capsettes™, Lime	500	10
6063	6043	Uni-Capsettes™/Capsettes™, Rose	500	10
6064	6044	Uni-Capsettes™/Capsettes™, Peach	500	10
6065	6045	Uni-Capsettes™/Capsettes™, Lavender	500	10



Capsettes™, Green 6033



Uni-Capsettes™, Yellow, 6051

**BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS**  
4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA  
Tel: (415) 472-3777 Fax: (415) 472-3758  
info@bioplas.com www.bioplas.com

*Histo Plas Capsettes™ and Uni-Capsettes™ are trademarks of Bio Plas, Inc., San Rafael, CA*