



BIO PLAS INC.

Premium Laboratory Consumables
Made In The USA



Pipet Tips: Reference Tips

Reference Pipet Tips				Reference Slim Pipet Tips			
A – 1-250 uL		B – 101-1000 uL		C – 1-200 uL		D – 101-1000 uL	
Replacement for: Eppendorf® 1-100 uL Pipetman™ 1-200 uL Finnpiptette 1-250 uL Socorex 5-100 uL		Replacement for: Oxford® 200-1000 uL Eppendorf® 101-1000 uL Pipetman™ 200-1000 uL Finnpiptette 200-1000 uL Socorex 101-1000 uL		Replacement for: Oxford® Series: 8000 10-200 uL 9000 Model Q 10-200 uL P-7000 1-200 uL 3000 Adjustable Abbott 1-200 uL Helena Quickpette 10-200 uL		Replacement for: Oxford® Series: 8000 600-1000 uL 9000 Model Q 600-1000 uL P-7000 200-1000 uL 3000 Adjustable Abbott 600-1000 uL Helena Quickpette 501-1000 uL	

Autoclavable
at
250°F
at
15 psi
for
15 mins

Ordering Information

A				B				C				D			
Catalog #	Tips/Rack	Tips/Pkg	Pkgs/Cs	Catalog #	Tips/Rack	Tips/Pkg	Pkgs/Cs	Catalog #	Tips/Rack	Tips/Pkg	Pkgs/Cs	Catalog #	Tips/Rack	Tips/Pkg	Pkgs/Cs
0001	N/A	1000	10	0002	N/A	1000	10	1012	N/A	1000	10	1212	N/A	1000	10
0001R	250	1000	10	0002R	250	1000	10	1012R	250	1000	10	1212R	250	1000	10
0001RS1	96	960	10	0002RS1	100	1000	5	N/A	N/A	N/A	N/A	1212RS1	100	1000	5

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com



Uni-Tip™ Universal Pipet Tips

Racked pipet tips allow the user to mount pipet tips in one efficient motion. The rack design features a seamless bottom that does not “breathe in” bench debris and contaminate tips when rack lid is opened. Racks are autoclavable at 250° F for 15 minutes at 15 psi.

1 – 250 uL		101 – 1000 uL	
	<p>Replacement for:</p> <p>Eppendorf® 1-100 uL Pipetman™ 1-200 uL Finnpiquette 1-250 uL Socorex 5-100 uL</p>		<p>Replacement for:</p> <p>Eppendorf® 101-1000 uL Pipetman™ 200-1000 uL Oxford® 200-1000 uL Finnpiquette 200-1000 uL Socorex 101-1000 uL</p>

Ordering Information							
1-250 uL				101-1000 uL			
Catalog #	Tips/ Rack	Tips/ Pkg	Pkg/Cs	Catalog #	Tips/ Rack	Tips/ Pkg	Pkg/Cs
0003	N/A	1000	10	0004	N/A	1000	10
0003R	250	1000	10	0004R	250	1000	10
0003RS1	96	960	10	0004S1	100	1000	10

- Uni-Tips™ Benefits:**
- Fit Multiple Brand Name Ejector Models
 - Inventory Reduction
 - Rimless
 - Racked or In Bulk

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 Info@bioplas.com www.bioplas.com



Economical Alternatives to Brand Name Pipet Tips

Made from virgin, non-wettable polypropylene, these tips are designed to replace specific manufacturer's tips.









For easy selection of replacement tips for your particular brand pipetter:

1. Find the brand you are currently using on the chart below. Example: Eppendorf®.
2. Find the tip volume. Example: 1-100 uL
3. Find the Tip Code at the bottom of the chart and use the code to locate the diagram and catalog number for the replacement tip in the charts on the following pages.

BRAND NAME CROSS REFERENCE																		
Manufacturer/Brand	Tip Volume Range, uL																	
C/A Selecapette	100- 1000																	
Eppendorf®		1-	101-															
Finnpipette				1- 250	200- 1000													
Gilson Pipetman™						1-	200-											
MLA®								10- 200	250- 1000									
Multi-Channel Titertek® and Octapette®											1- 250							
Oxford®											5- 200	1- 200	250- 500	600- 1000	250- 1000			
Oxford® Slim Tip												1- 200			250- 1000			
Sherwood/Lancer																5- 500		
Wheaton/Socorex																	5- 100	101- 1000
Tip Code	A	B	C	D	E	F	G	H	I	J	K	L	M	N	O	P	Q	R










BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com



REPLACEMENT TIP CHART			Catalog Numbers			
Tip Code	Diagram	Length (Inches)	Wells Per-Rack	Racks Per Pkg.	Racked Cat#	Bag of 1000 Cat#
A	 Natural	4	N/A	N/A	N/A	1800
B	 Yellow	2	250 96	4 10	1300R 1300RS1	1300
C	 Blue	3 1/2	250 100	4 10	1400R 1400RS1	1400
D	 Natural	2	250 96	4 10	2400R 2400RS1	2400
E	 Green	3 1/2	250 100	4 10	2500R 2500RS1	2500
F	 Yellow	2	250 96	4 10	2000R 2000RS1	2000
G	 Blue	3 1/2	250 100	4 10	2100R 2100RS1	2100
H	 Natural	1 1/2	250 96	4 10	1600R 1600RS1	1600

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com



REPLACEMENT TIP CHART			Catalog Numbers			
Tip Code	Diagram	Length (Inches)	Wells Per-Rack	Racks Per Pkg.	Racked	Bag of 1000
J		2	96	10	3500R	3500
			96	10	3500RS1	
K		1 ¼	250	4	1000R	1000
L		2 ½	250	4	1010R	1010
M		2 3/16				1100
N		3 ½	250	4	1200R	1200
			100	10	1200RS1	
O		3	250 100	4 10	1211R 1211RS1	1211
P		2 ¼				1500
Q		2	250	4	2600R	2600
			96	10	2600RS1	
R		1 ½	250 100	4 10	2700R 2700RS1	2700

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com



Ultra-Micro Pipet Tips

Designed to fit the full range of Eppendorf®, Pipetman™ and Oxford® Benchmate® Ultra Micro Pipettors. Thin, tapered delivery end eliminates carryover during pipetting. Manufactured with natural virgin polypropylene and warranted to be trace metal free. Racked pipet tips allow the user to mount pipet tips in one efficient motion. The rack design features a seamless bottom that does not “breathe in” bench debris and contaminate tips when rack lid is opened. Racks are autoclavable at 250° F for 15 minutes at 15 psi.

0.5 – 10 uL	0.1 – 10 uL
<div style="display: flex; align-items: center;"> <div> <p>Replacement for:</p> <p>Eppendorf® 0.5-10 uL</p> <p>Oxford® Benchmate® 0.5-10 uL</p> </div> </div>	<div style="display: flex; align-items: center;"> <div> <p>Replacement for:</p> <p>Pipetman™</p> <p>Model P-2 0.1-2.0 uL</p> <p>Model P-10 0.5-10 uL</p> </div> </div>


Ordering Information							
0.5 – 10 uL				0.1 – 10 uL			
Catalog #	Tips/Rack	Tips/Pkg	Pkg.Cs	Catalog #	Tips/Rack	Tips/Pkg	Pkg.Cs
1310	N/A	1000	10	2010	N/A	1000	10
1310RS1	96	960	10	2010RS1	96	960	10

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com



BIO PLAS INC.

Premium Laboratory Consumables
Made In The USA



Siliconized Pipet Tips

Universal fit that will provide an absolute seal on these single channel pipettors. *Siliconized Pipet Tips are not recommended for autoclaving.*

Siliconized Pipet Tips Benefits:

- Specially formatted siliconized polypropylene to help eliminate sample residue and cross-contamination
- Allows for precise and complete sample delivery
- Ideal for use in DNA and Thermal Cycling Procedures
- Low adhesion surface assures complete transfer of samples
- Prevents protein or nucleic acid adherence

1 – 200 uL		100 – 1000 uL	
	<p>Replacement for:</p> <p>Eppendorf® 1-100 uL Pipetman™ 1-200 uL Finnpiquette 1-250 uL Socorex 5-100 UI</p>		<p>Replacement for:</p> <p>Eppendorf® 101-1000 uL Pipetman™ 200-1000 uL Oxford® 250-1000 uL Finnpiquette 200-1000 uL Socorex 101-1000 uL</p>

Ordering Information

0.5 – 10 uL				0.1 – 10 uL			
Catalog #	Tips/Rack	Tips/Pkg	Pkg.Cs	Catalog #	Tips/Rack	Tips/Pkg	Pkg.Cs
3600SL	N/A	500	20	3605SL	N/A	500	20
3601SRS1	96	960	1	3606SRS1	100	1000	1

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com