



Bio Plas, Inc.

A Woman Owned Small Business

Premium Laboratory Consumables

Made In The USA

Uni-Flex™ Safety Caps Material Safety Data Sheet

Product Name: Uni-Flex™ Safety Caps

Bio Plas Catalog Numbers 6500 through 6730

Company Name: Bio Plas, Inc.

Address: 4380 Redwood Highway, Suite C-16, San Rafael, CA 94903

Tel: (415) 472-3777

Recommended Use: Capping Disposable Glass Culture Tubes and Blood Collection Tubes.

Product Composition: Low Density Polyethylene. Free of RNase, DNase, DNA, Pyrogens and PCR inhibitors.

Hazards Identification: Not classified as dangerous for use/supply

Hazardous Ingredients: None

Fire Hazard: Product burns at 349 Degrees Celsius. The fumes are CO and CO2

Melting Point: 110° C

Fire Extinguishing Media: Water, Foam, Carbon Dioxide, Dry Chemical

Toxicology Information: Material is non-toxic

Storage Requirements:

- This product may react with strong oxidising agents & should not be stored near such materials
- Storage temperature should be ambient. (preferably below 50 degrees C)
- Open flames prohibited.
- Protect from heat and direct sunlight.