



Bio Plas 5ml Flat Top Tubes

The perfect intermediate volume option between traditional microcentrifuge tubes and 15 mL conical tubes.

Bio Plas 5ml Flat Top Tube Features:

- Suitable for all rotors and racks that hold 15 mL conical tubes—no adapters needed.
- Graduated every 250 μ L for easy volume estimation.
- Rated to 25,000 x G for high-speed centrifugation.
- Flat, frosted, cap and a writing surface on the side of the tube.
- Temperature range: -80°C to 100°C.
- For use in pharmaceutical industry, food industry, and molecular biology.
- Manufactured from high-quality, virgin polypropylene to ensure purity and provide excellent optical clarity.
- Free of detectable RNase, DNase, DNA, PCR Inhibitors. Pyrogen-Free.
- Sterile tubes are certified endotoxin-free.



Cat# 4400



Cat# 4401



Cat# 4400AS

5ml Flat Top Tube Ordering Information

| Catalog # | Description | Color | Pack Qty | Packs/Case |
|-----------|---|----------|----------|------------|
| 4400 | 5ml Flat Top Microcentrifuge Tube | Clear | 200 | |
| 4401 | 5ml Flat Top Microcentrifuge Tube | Assorted | 200 | |
| 4400S | 5ml Flat Top Microcentrifuge Tube, Sterile, 4 Bags of 50 Per pack | Clear | 200 | |
| 4400AS | 5ml Flat Top Microcentrifuge Tube, Sterile, 4 Bags of 50 Per pack | Amber | 200 | |

Assorted colors include an assortment of 4 colors: blue, green, red, yellow. They come in 4 individually-packed, stand-up, resealable pouches for easy access and clean storage.

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16

San Rafael, CA 94903 USA

Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com www.bioplas.com