



BIO PLAS INC.
INNOVATIVE IDEAS...
FOR TODAY'S LABORATORY DISPOSABLES
MADE IN THE USA



Hexa-Flex™ Safety Caps

The Premium Multiple-Fit Tube Cap

A safe, economical way to cap and recap blood collection tubes and disposable glass culture tubes regardless of manufacturer.

Hexa-Flex™ Safety Caps are available for a wide range of biomedical and industrial applications and fit six different sized tubes. Hexa-Flex™ Safety Caps come in a variety of colors for coding samples.



The Hexa-Flex™ liquid-tight seal ensures sample integrity and prevents spillage of serum samples. Hexa-Flex™ Safety Caps are available in a variety of colors which helps guard against sample cross-contamination and are convenient for color-coding to particular R&D projects. Hexa-Flex™ Safety Caps protect your laboratory staff from the aerosols of highly-infectious organisms such as TB and the HTLV-III virus.



BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS

4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA

Tel: (415) 472-3777 Fax: (415) 472-3758

info@bioplas.com www.bioplas.com

Hexa-Flex™ is a trademarks of Bio Plas, Inc., San Rafael, CA



Hexa-Flex™ Safety Caps:

- Fit a variety of tubes, including 10,12, 13, 16 and 18mm sizes of blood collecting and disposable culture tubes.
- Liquid-tight seal prevents sample evaporation. The Hexa-Flex™ Safety Cap is designed to hold under centrifugation, agitation, refrigeration and freezing.
- Easy identification of samples. Hexa-Flex™ Safety Caps are available in seven colors that conveniently correspond to blood collection cap colors.
- Lock securely onto, or into, the tube lip to prevent evaporation.
- Easy, one-hand removal and application. To close, press down firmly with thumb to lock in place. To open, simply pry cap up.
- Clear, natural colored caps are trace-metal free.



Ordering Information			
Cat #	Description	Pkg. Qty	Pkgs./ Case
8350	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Natural	500	10
8355	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Red	500	10
8360	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Lavender	500	10
8365	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Green	500	10
8370	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Blue	500	10
8375	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Grey	500	10
8380	For 10mm, 12mm, 13mm, and 16mm and 18mm Blood Collection and Culture Tubes, Yellow	500	10

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com
 Hexa-Flex™ is a trademarks of Bio Plas, Inc., San Rafael, CA