



Bio Plas, Inc.

A Woman Owned Small Business

UNI-FLEX™ SAFETY CAPS

The Original Cap For Disposable Glass Culture And Blood Collection Tubes


A safe, economical, way to cap and recap blood collection tubes and disposable glass culture tubes regardless of manufacturer.

Uni-Flex™ Safety Caps are available for a wide range of biomedical and industrial applications.

Uni-Flex™ Safety Caps come in a variety of colors for coding samples.

The Uni-Flex™ liquid-tight seal ensures sample integrity and prevents spillage of serum samples.

Uni-Flex™ Safety Caps protect your laboratory staff from the aerosols of highly-infectious organisms such as TB and the HTLV-III virus.



Made In The USA!



Call Us
415-472-3777



Visit Our Website
www.bioplas.com



Bio Plas, Inc.

A Woman Owned Small Business

UNI-FLEX™ SAFETY CAPS

Features:

- 10, 12, 13 and 16mm sizes fit all varieties of blood collecting and disposable culture tubes.
- Will not aerosol when opened.
- Liquid-tight seal prevents sample evaporation. The Uni-Flex™ Safety Cap is designed to hold under centrifugation, agitation, refrigeration and freezing.
- Easy identification of samples. Uni-Flex™ Safety Caps are available in seven colors that conveniently correspond to blood collection cap colors.
- Lock securely onto the tube lip to prevent evaporation.
- Easy, one-hand removal and application. To close, press down firmly with thumb to lock in place. To open, simply pry cap up.
- Clear, natural colored caps are trace-metal free.



**Easy identification of samples:
Uni-Flex® Safety Caps are available in
seven colors that conveniently
correspond to blood collection cap colors**

Made In The USA!



Call Us
415-472-3777



Visit Our Website
www.bioplas.com



Bio Plas, Inc.

A Woman Owned Small Business

Uni-Flex™ Safety Cap Ordering Information

Catalog #	Description	Color	Pkg. Qty.	Pkgs./ Case
6500	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	White	1,000	10
6505	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Red	1,000	10
6510	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Lavender	1,000	10
6515	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Green	1,000	10
6520	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Blue	1,000	10
6525	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Grey	1,000	10
6530	Uni-Flex®Safety Caps for 10mm Blood Collecting & Culture Tubes	Yellow	1,000	10
6550	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	White	1,000	10
6555	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Red	1,000	10
6560	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Lavender	1,000	10
6565	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Green	1,000	10
6570	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Blue	1,000	10
6575	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Grey	1,000	10
6580	Uni-Flex®Safety Caps for 12mm Culture Tubes & 13mm Blood Collecting Tubes	Yellow	1,000	10
6600	Uni-Flex®Safety Caps for 13mm Culture Tubes	White	1,000	10
6605	Uni-Flex®Safety Caps for 13mm Culture Tubes	Red	1,000	10
6610	Uni-Flex®Safety Caps for 13mm Culture Tubes	Lavender	1,000	10
6615	Uni-Flex®Safety Caps for 13mm Culture Tubes	Green	1,000	10
6620	Uni-Flex®Safety Caps for 13mm Culture Tubes	Blue	1,000	10
6625	Uni-Flex®Safety Caps for 13mm Culture Tubes	Grey	1,000	10
6630	Uni-Flex®Safety Caps for 13mm Culture Tubes	Yellow	1,000	10
6700	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	White	1,000	10
6705	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Red	1,000	10
6710	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Lavender	1,000	10
6715	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Green	1,000	10
6720	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Blue	1,000	10
6725	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Grey	1,000	10
6730	Uni-Flex®Safety Caps for 16mm Blood Collecting & Culture Tubes	Yellow	1,000	10

Made In The USA!



Call Us
415-472-3777



Visit Our Website
www.bioplas.com