





## 80 Well Micro Tube Racks

The 80 Well Micro Tube Racks are designed with 80 wells in a 16 x 5 configuration with each well numbered for convenient sample identification. Rack dimensions are 9”L x 2 1/16”D with a well diameter of 10.5mm.

Available in 9 colors. This assortment of colors allows lab technologists to color code specific projects and allows for identification of hazardous experiments.



Ordering Information			
Catalog #	Description	Pkg. Qty	Pkgs./Case
0060	80 Well Micro Tube Rack, Natural	5	4
0061	80 Well Micro Tube Rack, Assorted*	5	4
0062	80 Well Micro Tube Rack, Blue	5	4
0063	80 Well Micro Tube Rack, Green	5	4
0064	80 Well Micro Tube Rack, Orange	5	4
0065	80 Well Micro Tube Rack, Pink	5	4
0066	80 Well Micro Tube Rack, Yellow	5	4
0067	80 Well Micro Tube Rack, Lavender	5	4
0068	80 Well Micro Tube Rack, Red	5	4
0069	80 Well Micro Tube Rack, Black	5	4



Autoclavable  
 at  
 250°F  
 at  
 15 psi  
 for  
 15 mins

**BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS**  
 4380 Redwood Hwy, Ste C-16, San Rafael, CA 94903 USA  
 Tel: (415) 472-3777 Fax: (415) 472-3758  
 info@bioplas.com www.bioplas.com