



96 and 64 Well “Reversible” Micro Tube Racks

Both the 96 and 64 Well “Reversible” Micro Tube Racks are designed so that one side of the rack accommodates 0.5ml microcentrifuge tubes. The reverse side holds 1.5ml to 2.0ml Flat Top Microcentrifuge Tubes and Screw Cap Microcentrifuge Tubes.


The clear polypropylene cover on the 96 Well “Reversible” Micro Tube Rack easily snaps on to keep it firmly in place during transportation and storage. A recessed area in the top of the cover allows for secure stacking of racks one on top of each other. The snap on cover is designed to accommodate the extra height of our Screw Cap Microcentrifuge Tubes. Two frosted writing areas for sample identification.

Racks are autoclavable: 15 minutes, 15lbs. @ 120° C and will withstand freezer storage to -90°C.

Alpha-numeric markings, “A” thru “H” vertically, 1 thru 8 horizontally.



64 Well Micro Tube Rack, Orange. Cat# 0084



96 Well Micro Tube Rack, Yellow. Cat# 0096

Ordering Information				
64 Well Catalog #	96 Well Catalog #	Description	Pkg. Qty	Pkgs./ Case
0080	0090	Micro Tube Rack, Natural	5	4
0081	0091	Micro Tube Rack, Assorted	5	4
0082	0092	Micro Tube Rack, Blue	5	4
0083	0093	Micro Tube Rack, Green	5	4
0084	0094	Micro Tube Rack, Orange	5	4
0085	0095	Micro Tube Rack, Pink	5	4
0086	0096	Micro Tube Rack, Yellow	5	4
0087	0097	Micro Tube Rack, Lavender	5	4
0089	0099	Micro Tube Rack, Black	5	4
0050	0052	Tube Rack Cover, Natural	5	4



Autoclavable
 at
 250°F
 at
 15 psi
 for
 15 mins

BIO PLAS, INC. A WOMAN OWNED SMALL BUSINESS
 4380 Redwood Hwy, Suite C-16, San Rafael, CA 94903 USA
 Tel: (415) 472-3777 Fax: (415) 472-3758
 info@bioplas.com www.bioplas.com